

ARAM User Guide


Version	1.0	Application date	23.04.2014
Code	REA-PUB-MOP-0001	Revision date	29.04.2015
Applied to		ARAM 1.2 beta	
Written by		Approved by	
Name	Alexandre KORNMANN	Abdelkrim BELHAOUA	
Function	C++ Developer	Scientist coordinator	

Document description

This document provides an installation and user guide of the library ARAM (Augmented Reality for Application on Medical field library).

Summary

1. ARAM User Guide	3
I. Prerequisite	3
II. Setting up	4
III. First application	5
IV. How does ARAM work?.....	10
V. Multi tracking	12
VI. ARAM API.....	16
2. Revision history.....	17

				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 2 sur 17

1. ARAM User Guide

I. Prerequisite

a. CMake

CMake is a cross-platform free software program for managing the build process of software using a compiler-independent method. It is designed to support directory hierarchies and applications that depend on multiple libraries, and for use in conjunction with native build environments such as make, Apple's Xcode, and Microsoft Visual Studio.

For ARAM library, **CMake** ≥ 2.8 is required.

b. OpenCV

OpenCV (Open Source Computer Vision) is a library of programming functions mainly aimed at real-time computer vision.

For ARAM library, **OpenCV** $\geq 2.4.8$ is required.

<p>Altran Research MEDIC@/alTRAN</p>				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 3 sur 17

II. Setting up

This guide was tested on Windows 7 Enterprise SP1 32bits and 64bits.

In this document, the default workspace is the directory **C:\ARAM**.

In this directory, we find (see ARAM WorkspaceFigure 1):

- **include**: a subdirectory for ARAM headers files (.hpp)
- **src**: a subdirectory for sources files (.cpp)
- **build**: an empty subdirectory, for builed files (.sln, .obj, ...)
- **CMakeLists.txt**: CMake configuration file, which help us to generate project file like **makefile** (gcc, clang, ...) or **Visual studio solution** (Visual studio).




Figure 1 ARAM Workspace

First, we generate our project file (Visual Studio solution, .sln), using **CMake GUI** (Figure 2).




Figure 2 CMake GUI configuration

If CMake fails in OpenCV linking, you have to indicate where OpenCV is in OpenCV_DIR variable.

After click on **Generate**, your directory should contain many files. Open **ARAM.sln**, and start the program!

<p>Altran Research MEDIC@/alTran</p>				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 4 sur 17

III. First application

In this section, our goal is to develop a simple application using ARAM library (Figure 3).




Figure 3 First application output

c. CMakeLists.txt

First, build your workspace for your application (as previously, **src**, **include**, **build**, and **CMakeLists.txt**.)

Typical, **CMakeLists.txt** looks as following:

```
cmake_minimum_required(VERSION 2.8)

project(ARAMProject)
set(EXECUTABLE_OUTPUT_PATH bin)

# ARAM include
include_directories("C:/ARAM/include")
link_directories("C:/ARAM/build/lib/Debug")

# OpenCV include
find_package(OpenCV REQUIRED)
include_directories(${OpenCV_INCLUDE_DIRS})

# Sources files
file(GLOB_RECURSE source_files src/* include/*)

add_executable(${CMAKE_PROJECT_NAME} ${source_files})

# linker
target_link_libraries(${CMAKE_PROJECT_NAME} ARAM)
target_link_libraries(${CMAKE_PROJECT_NAME} ${OpenCV_LIBS})
```

<p>Altran Research MEDIC@/alTran</p>				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 5 sur 17

d. Tags detection

This is a typical structure for ARAM application:

```
#include <ARAM/TagDetector.hpp> // Main ARAM class
#include <ARAM/tag/StandardTag.hpp> // Tag validator
#include <ARAM/ROIDetector/EdgeDetector.hpp> // Region of interest detection

#include <opencv2/opencv.hpp> // OpenCV data structure

#include <exception> // std::exception

int main(int argc, char** argv)
{
    try
    {
        // Detection parameters:
        // -> Region of interest detection
        // -> Tag validator
        typedef aram::TagDetector<aram::EdgeDetector, aram::StandardTag> myDetector;

        // Tag detector instantiation
        myDetector *detector = new myDetector();

        // Intrinsic parameters
        aram::Intrinsics intr("C:\\camera_data.xml");

        // Video input (see openCV doc)
        cv::VideoCapture cap(0); // use default video
        if(!cap.isOpened()) throw std::exception();

        cv::Mat frame;

        // Main loop
        while(true)
        {
            // next frame from video input
            cap >> frame;

            // Tag detection
            detector->detect(frame);

            // Tag list iterator
            aram::iteratorTag it;
            // Loop over valid tag in current frame
            for(it=detector->begin(); it!=detector->end(); ++it)
            {
                // Some operations here!
            }

            // render
            cv::imshow("render", frame);
            // GUI refresh (see openCV doc)
            if(cv::waitKey(10)>=0) break;
        }
    }
}
```

MEDIC@/ALTRAN				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 6 sur 17

```
catch(std::exception &)
{
}

return 0;
}
```

This program allows us to launch the camera.

Now, we add some operations using detected tag. First, we circle tags corners.
In this loop:

```
// Loop over valid tag in current frame
for(it=detector->begin();it!=detector->end();++it)
{
    // Some operations here!
}
```

As an example:

```
aram::vecPoint2D imgPoint = (*it)->corners();
for(unsigned int i=0;i<imgPoint.size();++i)
{
    cv::circle(frame,imgPoint[i],3,cv::Scalar(200,50,50),3);
}
```

Now, your output looks as following:




Figure 4 First application output (step 1)

Altran Research MEDIC@/alTran				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 7 sur 17

e. Pose estimation

For augmented reality, we need to compute the camera pose estimation. You can perform this with ARAM.

First, you have to calibrate your camera using OpenCV procedure. This procedure produces an .xml file, as the following:

I stored my .xml file at C:\camera_data.xml

```
<?xml version="1.0"?>
<opencv_storage>
<calibration_Time>"03/12/14 13:56:10"</calibration_Time>
<nrOfFrames>5</nrOfFrames>
<image_Width>1280</image_Width>
<image_Height>720</image_Height>
<board_Width>13</board_Width>
<board_Height>13</board_Height>
<square_Size>3.</square_Size>
<FixAspectRatio>1.</FixAspectRatio>
<!-- flags: +fix_aspectRatio +fix_principal_point +zero_tangent_dist -->
<flagValue>14</flagValue>
<Camera_Matrix type_id="opencv-matrix">
  <rows>3</rows>
  <cols>3</cols>
  <dt>d</dt>
  <data>
    1.2945600566643177e+003 0. 6.3950000000000000e+002 0.
    1.2945600566643177e+003 3.5950000000000000e+002 0. 0. 1.</data></Camera_Matrix>
<Distortion_Coefficients type_id="opencv-matrix">
  <rows>5</rows>
  <cols>1</cols>
  <dt>d</dt>
  <data>
    2.0026392183907174e-001 -2.7114152188654712e+000 0. 0.
    8.1223710401493587e+000</data></Distortion_Coefficients>
<Avg_Reprojection_Error>2.2073910408958071e+000</Avg_Reprojection_Error>
<Per_View_Reprojection_Errors type_id="opencv-matrix">
  <rows>5</rows>
  <cols>1</cols>
  <dt>f</dt>
  <data>
    1.11958861e+000 2.25522900e+000 3.74346900e+000 1.43035901e+000
    1.40137517e+000</data></Per_View_Reprojection_Errors>
</opencv_storage>
```

Based on previous code, add before the main loop this instruction:

```
// Intrinsics parameters
aram::Intrinsics intr("C:\\camera_data.xml");
```

And in the tag iteration loop:

Altran Research MEDIC@/ALTRAN				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 8 sur 17


```
// Loop over valid tag in current frame
for(it=detector->begin();it!=detector->end();++it)
{
    aram::vecPoint2D imgPoint = (*it)->corners();
    for(unsigned int i=0;i<imgPoint.size();++i)
    {
        cv::circle(frame,imgPoint[i],3,cv::Scalar(200,50,50),3);
    }

    // 3D points corresponding to corners
    aram::vecPoint3D objPoints;
    objPoints.push_back(aram::Point3D(0.0,0.0,0.0));
    objPoints.push_back(aram::Point3D(1.0,0.0,0.0));
    objPoints.push_back(aram::Point3D(1.0,1.0,0.0));
    objPoints.push_back(aram::Point3D(0.0,1.0,0.0));

    aram::Extrinsics e(intr,imgPoint,objPoints);

    aram::Point2D o = e.project(aram::Point3D(0.0,0.0,0.0));
    aram::Point2D x = e.project(aram::Point3D(1.0,0.0,0.0));
    aram::Point2D y = e.project(aram::Point3D(0.0,1.0,0.0));
    aram::Point2D z = e.project(aram::Point3D(0.0,0.0,1.0));

    cv::line(frame,o,x,cv::Scalar(255,0,0),3);
    cv::line(frame,o,y,cv::Scalar(0,255,0),3);
    cv::line(frame,o,z,cv::Scalar(0,0,255),3);
}
```

aram::Extrinsics compute the rotation matrix and translation vector between camera and tag, using openCV cv::solvePnP method. It allow you to project a 3D point (in world coordinate) to a 2D points (in image coordinate) using aram::Extrinsics::project method.

Be careful to the Point3D order!

Standard tag class order Point2D corners in "OpenCV order":

- Top left
- Bottom left
- Bottom right
- Top right





Figure 5 First application output (step 2)

Altran Research MEDIC@/ALTRAN				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 9 sur 17

IV. How does ARAM work?

f. Basics



Altran Research MEDIC@/alTRAN				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 10 sur 17

g. ARAM tag dictionary

ARAM is provided with a dictionary of 50 tags. tag[id].png, with id=0..49



Figure 6 tag0.png



Figure 7 tag1.png



Figure 8 tag3.png

<p>Altran Research MEDIC@/alTRAN</p>				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 11 sur 17

V. Multi tracking

ARAM can use a set of markers that define only one position! Multi tracking is a way to deal with occlusions.




Figure 9 Multi tracking output




Figure 10 Deal with occlusion

Altran Research MEDIC@/alTRAN				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 12 sur 17

h. Build a Grid

ARAM supports only coplanar tags detection.



Be consistent with your units! In this case, we use millimetres.

In your code:

```
float size = 142.0;
float delta = 654.0;

aram::Grid g;
aram::TagInfo t1(21,aram::Point2D(0.0,0.0),size);
aram::TagInfo t2(22, aram::Point2D(0.0,delta),size);
aram::TagInfo t3(19, aram::Point2D(delta,delta),size);
aram::TagInfo t4(20, aram::Point2D(delta,0.0),size);

g.addTagInfo(t1);
g.addTagInfo(t2);
g.addTagInfo(t3);
g.addTagInfo(t4);
```

Altran Research MEDIC@/alTran				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 13 sur 17

i. Multi tracking application

For multi tracking operations, we need to know tag's id. `aram::HammingTag` give an unique id, so we use `aram::HammingTag` instead of `aram::StandardTag` in our detector.

`aram::Chessboard::compute` return an `aram::Extrinsics`. We already use `aram::Extrinsics` in e. Pose estimation.

The minimal code for multi tracking is:

```
#include <ARAM/TagDetector.hpp> // Main ARAM class
#include <ARAM/tag/HammingTag.hpp> // Tag validator
#include <ARAM/ROIDetector/EdgeDetector.hpp> // Region of interest detection

#include <ARAM/tools/Grid.hpp> // Grid an TagInfo
#include <ARAM/coordinate/Chessboard.hpp> // Multi tracking Extrinsics compute

#include <opencv2/opencv.hpp> // OpenCV data structure

#include <exception> //std::exception

int main(int argc, char** argv)
{
    try
    {
        float size = 142.0;
        float delta = 654.0;

        aram::Grid g;
        aram::TagInfo t1(21,aram::Point2D(0.0,0.0),size);
        aram::TagInfo t2(22,aram::Point2D(0.0,delta),size);
        aram::TagInfo t3(19,aram::Point2D(delta,delta),size);
        aram::TagInfo t4(20,aram::Point2D(delta,0.0),size);

        g.addTagInfo(t1);
        g.addTagInfo(t2);
        g.addTagInfo(t3);
        g.addTagInfo(t4);
        aram::Chessboard *coord = new aram::Chessboard(g);

        // Detection parameters :
        // -> Region of interest detection
        // -> Tag validator USE HAMMINGTAG FOR MULTI TRACKING !
        typedef aram::TagDetector<aram::EdgeDetector,aram::HammingTag> myDetector;

        // Tag detector instantiation
        myDetector *detector = new myDetector();

        // Intrinsics parameters
        aram::Intrinsics intr("C:\\camera_data.xml");

        // Video input (see openCV doc)
        cv::VideoCapture cap(0); // use default video (usually your webcam)
```

```

if(!cap.isOpened()) throw std::exception();

cv::Mat frame;

// Main loop
while(true)
{
    // next frame from video input
    cap >> frame;

    // Tag detection
    detector->detect(frame);

    // Intrinsics parameters
    aram::Intrinsics intr("C:\\camera_data.xml");

    // Tag list iterator
    aram::iteratorTag it;

    // If any tags was detected
    if(detector->begin()!=detector->end())
    {
        // Get extrinsics parameters
        aram::Extrinsics e = coord->compute(detector->begin(),detector->end(),intr);

        // Project 3D world coordinate -> 2D image coordinate
        aram::Point2D o = e.project(aram::Point3D(0.0,0.0,0.0));
        aram::Point2D x = e.project(aram::Point3D(delta,0.0,0.0));
        aram::Point2D y = e.project(aram::Point3D(0.0,delta,0.0));
        aram::Point2D z = e.project(aram::Point3D(0.0,0.0,delta/2.0));

        // draw axis
        cv::line(frame,o,x,cv::Scalar(200,0,0),2);
        cv::line(frame,o,y,cv::Scalar(0,200,0),2);
        cv::line(frame,o,z,cv::Scalar(0,0,200),2);
    }

    // render
    cv::imshow("render", frame);
    // GUI refresh (see openCV doc)
    if(cv::waitKey(10)>=0) break;
}
}
catch(std::exception &)
{
}


return 0;
}

```

Altran Research MEDIC@/ALTRAN				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 15 sur 17

VI. ARAM API

j. UML diagram



2. Revision history

Revision history			
Version	Date	Author	Object
1.0	23.04.2014	Alexandre Kornmann	Creation

<p>Altran Research MEDIC@/alTRAN</p>				
Title:	ARAM User Guide		Code	REA-PUB-MOP-0001
Written by:	Alexandre Kornmann	Function:	C++ Developer	
Approved by:	Abdelkrim Belhaoua	Function:	Scientific coordinator	
Version :	1.0	Application date	23.04.2014	Page 17 sur 17